

**IMPORTANT SAFETY INFORMATION: Do not use the CryoGlow Mask if you are taking:**

- Medication or supplements that cause sensitivity to light (photosensitivity), for example, medications for severe inflammatory acne (i.e. isotretinoin or doxycycline).

**Immediately stop using the device if you experience any discomfort.**

**Do not use the CryoGlow Mask if you have any of the following medical conditions:**

- Retinal and/or macular disease, ocular albinism, or any eye condition that causes sensitivity to light
- Pregnant, nursing, or maybe pregnant
- Pacemaker or other electronic implanted device
- Implants in the face
- Lack of sensation (e.g., numbness)
- Light sensitive skin disorder (e.g., lupus, porphyria)
- Extreme light sensitivity of the eyes (i.e., photophobia)
- Light-induced headaches
- Albinism
- Melasma (i.e., condition that causes patches and spots)
- Sensitivity or allergy to plastic, aluminum, or self-fastening material

**Do not use the CryoGlow Mask if you have any of the following conditions on the face, hairline, ears, or neck:**

- Recent facial surgery (until you have been cleared by your doctor)
- Suspicious lesions or skin cancer
- Eczema
- Psoriasis
- Irregularly shaped, new, or changing brown spots
- Birthmarks
- Broken skin, skin rash, or blisters
- Tissue inflammation
- Skin infection, cold sore, or warts
- Active current outbreak of herpes simplex virus (cold sores)
- History of herpes outbreaks without receiving preventative treatment
- Facial treatment where the skin has not fully healed, including neurotoxin, dermal filler, micro needling, laser, or chemical peel

**Do not apply CryoGlow Mask's cooling routines on the eyeball or upper eyelid.**

**Do not use the CryoGlow Mask's cooling routines if you have any of the following medical conditions:**

- Cold hypersensitivity
- Cold urticaria (i.e., hives brought on by cold temperatures)
- Circulatory insufficiency

**PRECAUTIONS: Consult a medical professional before using the CryoGlow Mask if you:**

- Have bruising on the face
- Are experiencing skin sensitivity or irritation with retinol or benzoyl peroxide